

RAIN Sheath[®]

TIBIAL PEDAL INTRODUCER



Expanded Indications.
Expanded Access Options.

Cordis[®]

RAIN Sheath[®]

TIBIAL PEDAL INTRODUCER



HEXACUSPID HEMOSTASIS VALVE

Designed to preserve hemostasis and reduce risk of bleedback.



IN-VESSEL STABILITY

1 cm non-slip secure zone at proximal end of sheath designed to secure placement after insertion.




PROPRIETARY KINK RECOVERY TECHNOLOGY™

Elastometric properties allow the RAIN Sheath[®] Introducer to bend and flex to maintain lumen integrity.



LUBRICIOUS HYDROPHILIC COATING

Facilitates smoother, easier insertion and removal.



THIN-WALLED SHEATH

Available in 4F and 5F sheath sizes, designed for consistent performance during tibial pedal access.

RAIN Sheath[®]

TIBIAL PEDAL INTRODUCER

RAIN Sheath[®] Introducer for Tibial Pedal Procedures

SHEATH SIZE	CANNULA LENGTH	MINI-WIRE	MINI-WIRE COMPATABILITY	NEEDLE	PRODUCT CODE
4F	10 cm	Bare Stainless Steel	0.021" x 45 cm	21G Bare Needle	TP506410S
		Bare Nitinol	0.021" x 43 cm	21G Bare Needle	TP506410N
5F	10 cm	Bare Stainless Steel	0.021" x 45 cm	21G Bare Needle	TP506510S
		Bare Nitinol	0.021" x 43 cm	21G Bare Needle	TP506510N

Contraindications and Warnings for the RAIN Sheath[®] Tibial Pedal Introducer

INDICATIONS FOR USE

- The RAIN Sheath[®] Tibial Pedal Introducer is indicated to facilitate placing a catheter through the skin into the lower extremity peripheral vasculature below the knee.

CONTRAINDICATIONS

- None Known.

WARNINGS

- Use of alcohol, antiseptic solutions, or other solvents should be avoided, as they may adversely affect the device.
- Do not leave the CSI in place for extended periods of time without a catheter in place.
- Do not use a hydrophilic or polymer wire, with the bare needle, as this may damage the integrity of the wire coating or jacket.
- Manipulate the mini-guidewire slowly and carefully to avoid damage to the vessel wall, while monitoring the tip position and movement using standard catheterization technique.
- Once the vessel dilator is removed, manipulate the sheath introducer slowly and carefully to minimize the chances of kinking.
- Persons with allergic reactions to nickel may suffer an allergic response to components of this device.
- During the procedure, provide a proper anticoagulant or antiplatelet therapy to the patient.
- Do not use power injector for contrast media injection from the side port.
- Do not manually re-shape the tip of the mini-guidewire by applying external force intended to bend or affect the shape of mini-guidewire.
- This product is designed and intended for single use. It is not designed to undergo reprocessing and resterilization after initial use. Reuse of this product, including after reprocessing and/or resterilization, may cause a loss of structural integrity which could lead to a failure of the device to perform as intended and may lead to a difficult to clean after exposure to biological materials and may cause adverse patient reactions if reused.
- This device contains the following substance defined as CMR 1 A and/or CMR 1 B and/or endocrine disrupting substances in a concentration above 0.1 % weight by weight: Cobalt; CAS No. 7440-48-4; EC No. 231-158-0.
- Failure to abide by the warnings in this labeling might result in damage to the device coating, which may necessitate intervention or result in serious adverse events.

CAUTION: Federal (USA) law restricts this device to sale by or on the order of a physician.

Important information: Prior to use, refer to the instructions for use supplied with this device for indications, contraindications, side effects, suggested procedure, warnings and precautions.

CORDIS, Cordis LOGO, KINK RECOVERY TECHNOLOGY, and RAIN Sheath are trademarks of Cordis and may be registered in the US and/or in other countries.

©2023 Cordis. All Rights Reserved. 100630044 06/2023

Cordis[®]