



# Cordis S.M.A.R.T™

Nitinol Stent System

120 & 150  
mm length

## Think SMART



ALL LENGTHS COMPATIBLE WITH 6F SHEATH



# The longer S.M.A.R.T.™ Nitinol Stent System

Since their introduction in 1999, S.M.A.R.T.™ Stents have built and expanded their legacy. Every S.M.A.R.T.™ Stent is a reaffirmation of the Cordis Endovascular commitment to better patient care through evidence-based, results-oriented innovation.

**WE'RE NOT JUST INCREASING THE LENGTH.**


**WE'RE GROWING THE LEGACY.**

**120 & 150**  
mm length



## A legacy of proven design

- **EXCELLENT SCAFFOLDING**  
for uniform coverage
- **OFFSET ALTERNATING BRIDGES**  
for wall apposition on a bend
- **SHORT STRUT LENGTH**  
for individual cell response



## A legacy of confidence

- **MORE THAN 1,000,000 STENT PLACEMENTS WORLDWIDE**
- **ROBUST CLINICAL EVIDENCE**<sup>(1-2-3)</sup>
- **EXTENSIVELY STUDIED NITINOL STENT**<sup>(1-2-3)</sup>

<sup>1</sup>SIROCCO II, J. Endovasc. Ther. 2006; 13:701-710  
<sup>2</sup>CRISP\_US, J. Vasc. Interv. Radiol. 2004; 15:911-918  
<sup>3</sup>STROLL- 3 year results, 2013

Think **SMART**

# Cordis S.M.A.R.T.™ Nitinol Stent System

## Long 120 mm and 150 mm Stent

## Ordering Information

### Cordis S.M.A.R.T.™ Nitinol Stent System 120 mm and 150 mm Stent Lengths

120cm Delivery System	Expanded Stent Diameter (mm)	Expanded Stent Length (mm)	Profile: Sheath/Guide Fit (F)
C06120MV	6	120	6/8
C07120MV	7	120	6/8
C08120MV	8	120	6/8
C06150MV	6	150	6/8
C07150MV	7	150	6/8
C08150MV	8	150	6/8

### Cordis S.M.A.R.T. CONTROL™ Nitinol Stent System Large, 12 x 14 mm diameter

80cm Delivery System	120cm Delivery System	Expanded Stent Diameter (mm)	Expanded Stent Length (mm)	Profile: Sheath/Guide Fit (F)
C12030SV		12	30	7/9
C12040SV	C12040MV	12	40	7/9
C12060SV	C12060MV	12	60	7/9
C12080SV	C12080MV	12	80	7/9
C14030SV		14	30	7/9
C14040SV	C14040MV	14	40	7/9
C14060SV	C14060MV	14	60	7/9
C14080SV	C14080MV	14	80	7/9

### Cordis S.M.A.R.T. CONTROL™ Nitinol Stent System

80cm Delivery System	120cm Delivery System	Expanded Stent Diameter (mm)	Expanded Stent Length (mm)	Profile: Sheath/Guide Fit (F)
C06020SV		6	20	6/8
C06030SV	C06030MV	6	30	6/8
C06040SV	C06040MV	6	40	6/8
C06060SV	C06060MV	6	60	6/8
C06080SV	C06080MV	6	80	6/8
C06100SV	C06100MV	6	100	6/8
C07020SV		7	20	6/8
C07030SV	C07030MV	7	30	6/8
C07040SV	C07040MV	7	40	6/8
C07060SV	C07060MV	7	60	6/8
C07080SV	C07080MV	7	80	6/8
C07100SV	C07100MV	7	100	6/8
C08020SV		8	20	6/8
C08030SV	C08030MV	8	30	6/8
C08040SV	C08040MV	8	40	6/8
C08060SV	C08060MV	8	60	6/8
C08080SV	C08080MV	8	80	6/8
C08100SV	C08100MV	8	100	6/8
C09020SV		9	20	6/8
C09030SV	C09030MV	9	30	6/8
C09040SV	C09040MV	9	40	6/8
C09060SV	C09060MV	9	60	6/8
C09080SV	C09080MV	9	80	6/8
C10020SV		10	20	6/8
C10030SV	C10030MV	10	30	6/8
C10040SV	C10040MV	10	40	6/8
C10060SV	C10060MV	10	60	6/8
C10080SV	C10080MV	10	80	6/8