

# PALMAZ GENESIS™ Stent on OPTA™ Pro PTA Catheter (Large)

## .035" Over-the-wire System

PALMAZ GENESIS™ Stent on OPTA™ Pro (Large) PTA Catheter

- Recommended inflation pressure : 8 atm

## Ordering Information

Balloon diameter x length (mm x mm)	Unexpanded Stent Length (mm)	Expanded Stent Length (mm)	Sheath Compatibility (F)	Guiding Compatibility (F)	Catalog Number	
					Length of Delivery System 80 cm	Length of Delivery System 135 cm
5 x 40	39	39	6	8 (.088")	PG3950PPS	PG3950PPX
5 x 60	59	59	6	8 (.088")	PG5950PPS	PG5950PPX
6 x 30	29	27	6	8 (.088")	PG2960PPS	PG2960PPX
8 x 30	29	26	7	9 (.098")		PG2980PPX
6 x 40	39	38	6	8 (.088")	PG3960PPS	PG3960PPX
6 x 60	59	59	6	8 (.088")	PG5960PPS	PG5960PPX
6 x 80	79	79	6	8 (.088")	PG7960PPS	PG7960PPX
7 x 30	29	26	6	8 (.088")	PG2970PPS	PG2970PPX
7 x 40	39	36	6	8 (.088")		PG3970PPX
7 x 60	59	57	6.5	8 (.088")	PG5970PPS	PG5970PPX
7 x 80	79	78	6.5	8 (.088")	PG7970PPS	PG7970PPX
8 x 40	39	36	7	9 (.098")	PG3980PPS	PG3980PPX
8 x 60	59	56	7	9 (.098")	PG5980PPS	PG5980PPX
8 x 80	79	76	7	9 (.098")	PG7980PPS	PG7980PPX
9 x 20	19	17	7	9 (.098")	PG1990PPS	PG1990PPX
9 x 26	25	22	7	9 (.098")	PG2590PPS	PG2590PPX
9 x 30	29	25	7	9 (.098")	PG2990PPS	PG2990PPX
9 x 40	39	35	7	9 (.098")	PG3990PPS	PG3990PPX
9 x 60	59	53	7	9 (.098")	PG5990PPS	PG5990PPX
9 x 80	79	75	7	9 (.098")	PG7990PPS	PG7990PPX
10 x 20	19	16	7	9 (.098")	PG1910PPS	PG1910PPX
10 x 26	25	22	7	9 (.098")		PG2510PPX
10 x 30	29	24	7	9 (.098")	PG2910PPS	PG2910PPX
10 x 40	39	32	7	9 (.098")	PG3910PPS	PG3910PPX
10 x 60	59	53	7	9 (.098")	PG5910PPS	PG5910PPX
10 x 80	79	73	7	9 (.098")	PG7910PPS	

# Complementary Cordis Products

The PALMAZ GENESIS™ Stent on OPTA™ Pro PTA Catheter can be used with other Cordis products. The products listed hereunder are only a selection. Please contact your local sales representative for further information on product compatibilities.

## Ordering Information

### .035" STORQ™


Steerable Guidewires



Diameter (in)	Tip Flexibility	Tip Shape	Length (cm)	Catalogue Number
.035	Standard	Straight	180	503-456
.035	Standard	Angled	180	503-456J
.035	Standard	Modified J	180	503-456MJ
.035	Standard	Straight	300	503-456X
.035	Standard	Angled	300	503-456Y
.035	Standard	Modified J	300	503-456MY

### BRITE TIP™

Sheath Introducers



Cannula usable length (cm)	Catalogue Number 6F	Catalogue Number 7F	Catalogue Number 8F
11	401-611 M	401-711 M	401-811 M
23	401-623 M	401-723 M	401-823 M
35	401-635 M	401-735 M	401-835 M
45	401-645 M	401-745 M	401-845 M
55	401-655 M	401-755 M	401-855 M
90	401-690 M	401-790 M	401-890 M