

PALMAZ GENESIS™ Stent on OPTA™ Pro Catheter (Medium)

.035" Over-the-wire System

PALMAZ GENESIS™ Stent on OPTA™ Pro Catheter (Medium)

- Recommended inflation pressure: 8 atm

Ordering Information

Balloon diameter x length (mm x mm)	Unexpanded Stent Length (mm)	Expanded Stent Length (mm)	Sheath Compatibility (F)	Guiding Compatibility 8F / .088"	Catalog Number	
					Length of Delivery System 80 cm	135 cm
4 x 17	15	14	6	8	PG1540PPS	
4 x 20	18	17	6	8	PG1840PPS	
5 x 15	12	11	6	8	PG1250PPS	PG1250PPX
5 x 17	15	14	6	8	PG1550PPS	
5 x 17	15	14	6	8	PG1550PPS	
5 x 20	18	17	6	8	PG1850PPS	PG1850PPX
5 x 25	24	23	6	8	PG2450PPS	
6 x 17	15	14	6	8	PG1560PPS	
6 x 20	18	17	6	8	PG1860PPS	PG1860PPX
6 x 25	24	23	6	8	PG2460PPS	PG2460PPX
7 x 15	12	11	6	8	PG1270PPS	PG1270PPX
7 x 17	15	14	6	8	PG1570PPS	
7 x 20	18	17	6	8	PG1870PPS	PG1870PPX
7 x 25	24	23	6	8	PG2470PPS	PG2470PPX
8 x 15	12	11	6	8	PG1280PPS	PG1280PPX
8 x 17	15	14	6	8	PG1580PPS	PG1580PPX
8 x 20	18	17	6	8	PG1880PPS	PG1880PPX
8 x 25	24	23	6	8	PG2480PPX	

Complementary Cordis Products

The PALMAZ GENESIS™ Stent on OPTA™ Pro Catheter can be used with other Cordis Products or Product Families. The product listed below is only a selection. please contact your local sales representative for further information on product compatibilities.

Ordering Information

VISTA BRITE TIP™ IG

Introducer Guides



Catheter Shape	Length (cm)	OD (F)	ID	Catalogue Number
Renal curve (small)	55	8	.088"	403-8553R
Renal curve (long)	55	8	.088"	403-8553A

For Healthcare Professionals Only.

Important information: Prior to use, refer to the instructions for use supplied with this device for indications, contraindications, side effects, suggested procedure, warnings and precautions. As part of its continuous product development policy, Cordis reserves the right to change product specifications without prior notification. Please contact your Cordis representative for additional product availability information. CORDIS, Cordis LOGO, JINDO, VISTA BRITE TIP, OPTA and PALMAZ GENESIS are trademarks of Cordis and may be registered in the US and/or in other countries.

© 2022 Cordis. All Rights Reserved. 100549936-3 12/2022 .



Cordis™