

**Product Description**

**Type:** Low profile, high pressure, over-the-wire PTA balloon catheter designed for use in the iliac, femoral, ilio-femoral, popliteal, infra popliteal and renal arteries. The device is also indicated for post-dilatation of balloon expandable and self-expanding stents in the peripheral vasculature.

<b>RBP RANGE</b>	12-24 atm	<b>SHEATH FIT RANGE</b>	5F-7F
<b>NOMINAL PRESSURE</b>	8-10 atm	<b>BALLOON MATERIAL</b>	Nylon (Grilamid)
<b>BALLOON LENGTHS</b>	20-300 mm	<b>MARKERS</b>	Proximal and distal marker bands (dual distal marker bands for balloon lengths ≥100 mm)
<b>BALLOON DIAMETERS</b>	3-12 mm	<b>GUIDEWIRE</b>	.035"
<b>SHAFT LENGTHS</b>	80 cm, 135 cm	<b>COATING</b>	Dual layer hydrophilic coating
<b>SHAFT MATERIAL</b>	Pebax		



**Ordering Information** SABER™ .035 PTA DILATATION CATHETER

DIAMETER (mm)	LENGTH (mm)	RBP	SHEATH	SHAFT LENGTH	
				80 cm	135 cm
3	20	24	5	48003002N	48003002R
3	40	24	5	48003004N	48003004R
3	60	24	5	48003006N	48003006R
3	80	22	5	48003008N	48003008R
3	100	22	5	48003010N	48003010R
3	120	22	5	48003012N	48003012R
3	150	22	5	48003015N	48003015R
3	200	22	5	48003020N	48003020R
3	250	20	5	48003025N	48003025R
3	300	20	5	48003030N	48003030R

DIAMETER (mm)	LENGTH (mm)	RBP	SHEATH	SHAFT LENGTH	
				80 cm	135 cm
4	20	22	5	48004002N	48004002R
4	40	22	5	48004004N	48004004R
4	60	22	5	48004006N	48004006R
4	80	22	5	48004008N	48004008R
4	100	22	5	48004010N	48004010R
4	120	22	5	48004012N	48004012R
4	150	22	5	48004015N	48004015R
4	200	22	5	48004020N	48004020R
4	250	20	5	48004025N	48004025R
4	300	20	5	48004030N	48004030R

# Ordering Information

## SABER™ .035 PTA DILATATION CATHETER

DIAMETER (mm)	LENGTH (mm)	RBP	SHEATH	SHAFT LENGTH	
				80 cm	135 cm
5	20	22	5	48005002N	48005002R
5	40	22	5	48005004N	48005004R
5	60	22	5	48005006N	48005006R
5	80	22	5	48005008N	48005008R
5	100	22	5	48005010N	48005010R
5	120	20	5	48005012N	48005012R
5	150	20	5	48005015N	48005015R
5	200	20	5	48005020N	48005020R
5	250	18	5	48005025N	48005025R
5	300	18	5	48005030N	48005030R
6	20	18	5	48006002N	48006002R
6	40	18	5	48006004N	48006004R
6	60	17	5	48006006N	48006006R
6	80	17	5	48006008N	48006008R
6	100	17	5	48006010N	48006010R
6	120	17	5	48006012N	48006012R
6	150	17	5	48006015N	48006015R
6	200	17	6	48006020N	48006020R
6	250	17	6	48006025N	48006025R
6	300	17	6	48006030N	48006030R
7	20	18	5	48007002N	48007002R
7	40	18	5	48007004N	48007004R
7	60	16	5	48007006N	48007006R
7	80	16	5	48007008N	48007008R
7	100	16	6	48007010N	48007010R

DIAMETER (mm)	LENGTH (mm)	RBP	SHEATH	SHAFT LENGTH	
				80 cm	135 cm
7	120	14	6	48007012N	48007012R
7	150	14	6		48007015R
7	200	14	6		48007020R
7	250	14	6		48007025R
7	300	14	6		48007030R
8	20	16	6	48008002N	48008002R
8	40	16	6	48008004N	48008004R
8	60	16	6	48008006N	48008006R
8	80	16	6	48008008N	48008008R
8	100	16	6	48008010N	48008010R
9	20	13	6	48009002N	48009002R
9	40	13	6	48009004N	48009004R
9	60	13	6	48009006N	48009006R
9	80	13	6	48009008N	48009008R
9	100	13	6	48009010N	48009010R
10	20	14	7	48010002N	48010002R
10	40	14	7	48010004N	48010004R
10	60	14	7	48010006N	48010006R
10	80	14	7	48010008N	48010008R
10	100	14	7	48010010N	48010010R
12	20	12	7	48012002N	48012002R
12	40	12	7	48012004N	48012004R

For Healthcare Professionals Only.

**Important information:** Prior to use, refer to the instructions for use supplied with this device for indications, contraindications, side effects, suggested procedure, warnings and precautions. As part of its continuous product development policy, Cordis reserves the right to change product specifications without prior notification. Please contact your Cordis representative for additional product availability information.

CORDIS, Cordis LOGO, and SABER are trademarks of Cordis and may be registered in the US and/or in other countries.

© 2022 Cordis. All Rights Reserved. 100573563-2 08/2022