


# Cordis STABILIZER™ Support Guidewire and Cordis STABILIZER™ Extra Support Guidewire

0.014" Peripheral Steerable Guidewires



Providing All the Support  
You Need for 0.014"  
Peripheral Interventions

# Cordis STABILIZER™ Support Guidewire and Cordis STABILIZER™ Extra Support Guidewire

## 0.014" Peripheral Steerable Guidewires

### Versatility

Broad transitions and high support for device deliverability and compatibility with OUTBACK™ Elite Re-Entry Catheter

### Enhanced Lubricity

Proximal DURAGLIDE™ PTFE Spray and distal PTFE Sleeve for excellent surface lubricity

### Atraumatic Tip

SuperSoft tip flexibility provides atraumatic vessel entry and navigation

## Ordering Information

### Cordis STABILIZER™ Support Guidewire

Catalog Numbers	Diameter (inches)	Tip Flexibility	Length (cm)	Tip Shape
507 - 180S	0.014	SuperSoft	180	Straight
507 - 300S	0.014	SuperSoft	300	Straight

### Cordis STABILIZER™ Extra Support Guidewire


Catalog Numbers	Diameter (inches)	Tip Flexibility	Length (cm)	Tip Shape
527 - 180E	0.014	SuperSoft	180	Straight
527 - 300E	0.014	SuperSoft	300	Straight

### Compatible with OUTBACK™ Elite Re-Entry Catheter

Codes	Length (cm)	Sheath (F)	-	-
OTB 59080A	80	6	-	-
OTB 59120A	120	6	-	-

### Cordis STABILIZER™ Support Guidewire

Support for the delivery of advanced peripheral devices in complex and tortuous anatomy



### Cordis STABILIZER™ Extra-Support Guidewire

When a higher level of support is needed for the most challenging peripheral interventions

