

SUPER TORQUE™ MB

5F Marker Band Flush Catheters



Cordis™


SUPER TORQUE™ MB

5F Marker Band Flush Catheters

Matching the exacting standards of Cordis quality, the SUPER TORQUE™ MB markerband catheter was designed to provide angiographic visualization and linear measurement of the vasculature when combined with the delivery of radiopaque contrast media to selected sites in the vascular system.

The excellent visibility of the gold alloy marker bands provide easy location information of the catheter in the vasculature.

The braided body provides you with a manoeuvrable device that is demanded by this type of challenging procedure.



Ordering Information

- Sold by 5 units

Lengths (cm)	Size (F)	2 gold-alloy marker bands		10 gold-alloy marker bands		20 gold-alloy marker bands	
		SH*	Cat. N°	SH*	Cat. N°	SH*	Cat. N°
65	5	-		-		8	532-598C
110	5	-		6	532-598A	-	532-598B
125	5	6		-		-	SRD7040MB

*Side Holes