

Cordis[®]












BRITE TIP[®]

Guiding Catheter

ブライチップ[®]
ガイディングカテーテル

Reliable BLUE



先端形状		カタログ番号 [外径 (F) / 内径 (inch)]					
		有効長 (cm)	6F/0.070"	7F/0.078"	8F/0.088"	9F/0.098"	10F/0.110"
	Judkins Right4 (JR4)	55	670-082-55	—	—	—	—
		110	SM7504	—	—	—	—
	Internal Mammary (IM)	55	670-190-55	—	—	—	—
	Hockey Stick (HS)	55	670-278-55	—	588-841P*	—	—
	Multipurpose (MPA I)	90	670-270-90	778-270-90	588-842P	598-942P*	—
		120	SM7472	—	—	—	—
	Straight (STR)	55	670-280-55	—	—	—	—
		80	—	—	—	—	511-040P
		90	670-280-90	778-280-90	588-844P	598-944P*	—
	Renal Curve (Small) (RDC I)	55	670-210-55	—	588-845P*	—	—
	Renal Curve (RDC)	55	670-212-55*	—	588-846P*	—	—
	Vertebral (VTBL)	90	—	778-216-90	—	—	—
	J-Curve (J)	55	670-214-55*	—	—	—	—
	Multipurpose (MPD)	95	—	—	588-8207*	—	—
	Headhunter 1 (H1)	95	—	—	588-8202*	—	—

*印のついている製品は在庫をご確認の上、ご注文ください。

[7,8,9F]

■販売名：コーディス プライットチップ ガイディングカテーテル ■承認番号：20900BZY00298000 ■クラス分類：IV ■JMDNコード：17846104

■特定保険医療材料請求区分：ガイディングカテーテル/その他血管用

[6,10F]

■販売名：プライットチップ ガイディングカテーテル ■承認番号：21200BZY00636000 ■クラス分類：IV ■JMDNコード：17846104

■特定保険医療材料請求区分：ガイディングカテーテル/その他血管用

