


Product Description

Type: Low profile, high pressure, over-the-wire PTA balloon catheter designed for use in the iliac, femoral, ilio-femoral, popliteal, infra popliteal and renal arteries and for the treatment of obstructive lesions of native or synthetic arteriovenous dialysis fistulae. The device is also indicated for post-dilatation of balloon expandable and self-expanding stents in the peripheral vasculature.

RBP RANGE	12-24 atm	SHEATH FIT RANGE	5F-7F
NOMINAL PRESSURE	8-10 atm	BALLOON MATERIAL	Nylon (Grilamid)
BALLOON LENGTHS	20-300 mm	MARKERS	Proximal and distal marker bands (dual distal marker bands for balloon lengths ≥100 mm)
BALLOON DIAMETERS	3-12 mm	GUIDEWIRE	.035"
SHAFT LENGTHS	40 cm, 80 cm, 135 cm	COATING	Dual layer hydrophilic coating
SHAFT MATERIAL	Pebax		



Ordering Information **SABER® .035 PTA DILATATION CATHETER**

DIAMETER (mm)	LENGTH (mm)	RBP	SHEATH	SHAFT LENGTH		
				40 cm	80 cm	135 cm
3	20	24	5	48003002P	48003002N	48003002R
3	40	24	5	48003004P	48003004N	48003004R
3	60	24	5	48003006P	48003006N	48003006R
3	80	22	5	48003008P	48003008N	48003008R
3	100	22	5	48003010P	48003010N	48003010R
3	120	22	5		48003012N	48003012R
3	150	22	5		48003015N	48003015R
3	200	22	5		48003020N	48003020R
3	250	20	5		48003025N	48003025R
3	300	20	5		48003030N	48003030R

DIAMETER (mm)	LENGTH (mm)	RBP	SHEATH	SHAFT LENGTH		
				40 cm	80 cm	135 cm
4	20	22	5	48004002P	48004002N	48004002R
4	40	22	5	48004004P	48004004N	48004004R
4	60	22	5	48004006P	48004006N	48004006R
4	80	22	5	48004008P	48004008N	48004008R
4	100	22	5	48004010P	48004010N	48004010R
4	120	22	5		48004012N	48004012R
4	150	22	5		48004015N	48004015R
4	200	22	5		48004020N	48004020R
4	250	20	5		48004025N	48004025R
4	300	20	5		48004030N	48004030R

Ordering Information

SABER® .035 PTA DILATATION CATHETER

DIAMETER (mm)	LENGTH (mm)	RBP	SHEATH	SHAFT LENGTH		
				40 cm	80 cm	135 cm
5	20	22	5	48005002P	48005002N	48005002R
5	40	22	5	48005004P	48005004N	48005004R
5	60	22	5	48005006P	48005006N	48005006R
5	80	22	5	48005008P	48005008N	48005008R
5	100	22	5	48005010P	48005010N	48005010R
5	120	20	5	48005012P	48005012N	48005012R
5	150	20	5		48005015N	48005015R
5	200	20	5		48005020N	48005020R
5	250	18	5		48005025N	48005025R
5	300	18	5		48005030N	48005030R
6	20	18	5	48006002P	48006002N	48006002R
6	40	18	5	48006004P	48006004N	48006004R
6	60	17	5	48006006P	48006006N	48006006R
6	80	17	5	48006008P	48006008N	48006008R
6	100	17	5	48006010P	48006010N	48006010R
6	120	17	5	48006012P	48006012N	48006012R
6	150	17	5		48006015N	48006015R
6	200	17	6		48006020N	48006020R
6	250	17	6		48006025N	48006025R
6	300	17	6		48006030N	48006030R
7	20	18	5	48007002P	48007002N	48007002R
7	40	18	5	48007004P	48007004N	48007004R
7	60	16	5	48007006P	48007006N	48007006R
7	80	16	5	48007008P	48007008N	48007008R
7	100	16	6	48007010P	48007010N	48007010R

DIAMETER (mm)	LENGTH (mm)	RBP	SHEATH	SHAFT LENGTH		
				40 cm	80 cm	135 cm
7	120	14	6		48007012N	48007012R
7	150	14	6			48007015R
7	200	14	6			48007020R
7	250	14	6			48007025R
7	300	14	6			48007030R
8	20	16	6	48008002P	48008002N	48008002R
8	40	16	6	48008004P	48008004N	48008004R
8	60	16	6	48008006P	48008006N	48008006R
8	80	16	6	48008008P	48008008N	48008008R
8	100	16	6	48008010P	48008010N	48008010R
9	20	13	6	48009002P	48009002N	48009002R
9	40	13	6	48009004P	48009004N	48009004R
9	60	13	6	48009006P	48009006N	48009006R
9	80	13	6	48009008P	48009008N	48009008R
9	100	13	6	48009010P	48009010N	48009010R
10	20	14	7	48010002P	48010002N	48010002R
10	40	14	7	48010004P	48010004N	48010004R
10	60	14	7	48010006P	48010006N	48010006R
10	80	14	7	48010008P	48010008N	48010008R
10	100	14	7	48010010P	48010010N	48010010R
12	20	12	7	48012002P	48012002N	48012002R
12	40	12	7	48012004P	48012004N	48012004R