



Closure starts with a secure fit. The ZEPHYR® Vascular Compression Band is designed to provide security, comfort, and simplicity. This band is uniquely universal with its standard syringe compatibility and multiple access site versatility.



**SOFT FLEXIBLE BAND**

Compliant elastomeric band for patient comfort and firm compression balloon for patent hemostasis



**INTERCHANGABLE SYRINGE CONNECTION**

Easy to use: universally compatible with standard luer syringes



**DOUBLE BONDED RADIAL BALLOON**

The ZEPHYR® Band help clinicians achieve patent hemostasis with firm downward pressure and clear visualization of the puncture site


PRODUCT CODE	ZEPHYR BAND TYPE	LENGTH	UNITS
190102	Large	30cm	5EA/BOX

Part of the Radial360 portfolio, ZEPHYR® Band also works for multiple access points.

The ZEPHYR® Vascular Compression Band is indicated for use by medical professionals to promote hemostasis following a catheterization or other puncture into a blood vessel in a patient's arm or leg, including: radial, brachial, dorsalis pedis, or tibial blood vessels, arterial or venous line or sheath removal, hemodialysis, and in anticoagulation therapy.

CAUTION: Federal (USA) law restricts this device to sale by or on the order of a physician. Important information: Prior to use, refer to the instructions for use supplied with this device for indications, contraindications, side effects, suggested procedure, warnings and precautions. ZEPHYR® Vascular Compression Band is manufactured by Advanced Vascular Dynamics and distributed by Cordis. CORDIS and Cordis LOGO are trademarks of Cordis and may be registered in the US and/or in other countries. All other marks are the property of their respective owners. ©2021 Cordis. All Rights Reserved. 100525155-2 11/2021

# ONE BAND



Simplify Patent Hemostasis