


Cordis

BRITE TIP[®]

Sheath Introducers

—

Illuminating the Path to
Precision Placement



Kink resistant cannula

RESILIENT DUAL LAYER CANNULA


- Flexibility and support to access tortuous vessels

REINFORCED CANNULA WALL

- Kink resistance

SILICONE COATED CANNULA AND VESSEL DILATOR

- Lubricious, anti-thrombogenic vessel entry



SLIX valve

SIX-CUT SPIRAL DESIGN FOR HEMOSTATIC VALVE

- Free catheter movement without blood reflux or air aspiration




Manicured tip and dilator design

Cordis

TAPERED TRANSITIONS ON SHEATH AND DILATOR FOR:

- Atraumatic vessel entry
- Improved patient comfort



Radiopaque BRITE TIP[®] Sheath

BLACK TUNGSTEN-FILLED DISTAL TIP


- Excellent visualization facilitates exact location of sheath's distal tip for accurate interventional positioning

BRITE TIP®

Sheath Introducers

ORDERING INFORMATION

Longer lengths available
for a wide variety of
interventional procedures



CATALOG NUMBERS

Cannula Length			4.5cm	5.5cm	11cm	23cm	35cm	45cm	55cm	70cm	90cm
French Size	Color Code	Dilator Tip (mm)	With Mini-Wire	With Mini-Wire	With Mini-Wire	Without Mini-Wire	Without Mini-Wire	Without Mini-Wire	Without Mini-Wire	Without Mini-Wire	Without Mini-Wire
4	■	25	SRA0116		401-411M	401-423M	SRA0207	SRAM17661			
5	■	25	SRA0117	SRA0035	401-511M	401-523M	401-535M	401-545M	401-555M		
6	■	35	SRA0115	401-605M	401-611M	401-623M	401-635M	401-645M	401-655M	SRA0217	401-690M
7	■	35		401-705M	401-711M	401-723M	401-735M	401-745M	401-755M		401-790M
8	■	35		401-805M	401-811M	401-823M	401-835M	401-845M	401-855M		
9	■	35			401-911M	401-923M	401-935M				
10	■	45			401-011M	401-023M					
11	■	45			401-111M	401-123M					

*Each catalog number must be ordered by the box

QUANTITY

- 5 units per package for 4.5cm, 5.5cm, 11cm and 23cm lengths
- 2 units per package for 35cm or longer lengths

CAUTION: Federal (USA) law restricts this device to sale by or on the order of a physician.

Important information: Prior to use, refer to the instruction for use supplied with this device for indications, contraindications, side effects, suggested procedure, warnings and precautions.

CORDIS, Cordis LOGO, and BRITE TIP are trademarks of Cordis and may be registered in the US and/or in other countries.
© 2021 Cordis. All Rights Reserved. 100576607-2 12/2021

Cordis