

## PRODUCT ORDER FORM

### BRITE TIP RADIANTZ™ Guiding Sheath

The first multi-segmented thin-walled guiding sheath specifically designed to optimize lower-extremity lesion access from a radial approach.

Ordering Information (Packaging: 1 unit/box)



SIZE (Fr)	DILATOR ID	SHAFT LENGTH (cm)	SKU	ORDER QTY
6	.035	110	687280110	
6	.035	135	687280135	

### S.M.A.R.T. RADIANTZ™ Vascular Stent System

Purposely engineered to treat iliac, SFA and proximal popliteal lesions from the radial approach. The first radial peripheral self-expanding stent indicated for iliac lesions.

Ordering Information (Packaging: 1 unit/box)



UNCONSTRAINED STENT DIAMETER (mm)	UNCONSTRAINED STENT LENGTH (mm)	RECOMMENDED VESSEL DIAMETER (mm)	SHEATH COMPATIBILITY (Fr)	STENT INDICATION		SHAFT LENGTH		ORDER QTY
						Compatible with 110cm BRITE RADIANTZ™ Guiding Sheath	Compatible with 135cm BRITE TIP RADIANTZ™ Guiding Sheath	
				ILIAC	SFA	150cm	190cm	
6	20	4 - 5	6		x		SR06020XL	
6	30	4 - 5	6		x		SR06030XL	
6	40	4 - 5	6		x		SR06040XL	
6	60	4 - 5	6		x		SR06060XL	
6	80	4 - 5	6		x		SR06080XL	
6	100	4 - 5	6		x		SR06100XL	
6	120	4 - 5	6		x		SR06120XL	
6	150	4 - 5	6		x		SR06150XL	
7	20	5 - 6	6		x		SR07020XL	
7	30	5 - 6	6		x		SR07030XL	
7	40	5 - 6	6		x		SR07040XL	
7	60	5 - 6	6		x		SR07060XL	
7	80	5 - 6	6		x		SR07080XL	
7	100	5 - 6	6		x		SR07100XL	
7	120	5 - 6	6		x		SR07120XL	
7	150	5 - 6	6		x		SR07150XL	
8	20	6 - 7	6	x	x	SR08020L	SR08020XL	
8	30	6 - 7	6	x	x	SR08030L	SR08030XL	
8	40	6 - 7	6	x	x	SR08040L	SR08040XL	
8	60	6 - 7	6	x	x	SR08060L	SR08060XL	
8	80	6 - 7	6	x	x	SR08080L	SR08080XL	
8	100	6 - 7	6	x	x	SR08100L	SR08100XL	
8	120	6 - 7	6	x	x	SR08120L	SR08120XL	
8	150	6 - 7	6	x	x	SR08150L	SR08150XL	
9	20	7 - 8	6	x		SR09020L		
9	30	7 - 8	6	x		SR09030L		
9	40	7 - 8	6	x		SR09040L		
9	60	7 - 8	6	x		SR09060L		
9	80	7 - 8	6	x		SR09080L		
10	20	8 - 9	6	x		SR10020L		
10	30	8 - 9	6	x		SR10030L		
10	40	8 - 9	6	x		SR10040L		
10	60	8 - 9	6	x		SR10060L		
10	80	8 - 9	6	x		SR10080L		

## PRODUCT ORDER FORM

### SABERX RADIANZ™ Percutaneous Transluminal Angioplasty (PTA) Dilatation Catheter

The most clinically comprehensive radial peripheral PTA size matrix that facilitates broader treatment of lesions and delivers exceptional performance.

Ordering Information (Packaging: 1 unit/box)



BALLOON DIAMETER (mm)	BALLOON LENGTH (mm)	RATED BURST PRESSURE (atm)	SHEATH COMPATIBILITY (Fr)	SHAFT LENGTH (190cm)	ORDER QTY
2	20	18	4	R51002002L	
2	30	18	4	R51002003L	
2	40	18	4	R51002004L	
2	60	18	4	R51002006L	
2	80	18	4	R51002008L	
2	100	18	4	R51002010L	
2	150	18	4	R51002015L	
2	200	18	4	R51002020L	
2	250	18	4	R51002025L	
2	300	18	4	R51002030L	
2.5	20	18	4	R51002502L	
2.5	30	18	4	R51002503L	
2.5	40	18	4	R51002504L	
2.5	60	18	4	R51002506L	
2.5	80	18	4	R51002508L	
2.5	100	18	4	R51002510L	
2.5	150	18	4	R51002515L	
2.5	200	18	4	R51002520L	
2.5	250	18	4	R51002525L	
2.5	300	18	4	R51002530L	
3	20	18	4	R51003002L	
3	30	18	4	R51003003L	
3	40	18	4	R51003004L	
3	60	18	4	R51003006L	
3	80	18	4	R51003008L	
3	100	18	4	R51003010L	
3	150	14	4	R51003015L	
3	200	14	4	R51003020L	
3	250	14	4	R51003025L	
3	300	14	4	R51003030L	
3.5	20	18	4	R51003502L	
3.5	30	18	4	R51003503L	
3.5	40	18	4	R51003504L	
3.5	60	18	4	R51003506L	
3.5	80	18	4	R51003508L	
3.5	100	18	4	R51003510L	
3.5	150	14	4	R51003515L	
3.5	200	14	4	R51003520L	
3.5	250	14	4	R51003525L	
3.5	300	14	4	R51003530L	
4	20	18	4	R51004002L	
4	30	18	4	R51004003L	
4	40	18	4	R51004004L	
4	60	18	4	R51004006L	
4	80	18	4	R51004008L	

# PRODUCT ORDER FORM

## SABERX RADIANZ™ Percutaneous Transluminal Angioplasty (PTA) Dilatation Catheter

The most clinically comprehensive radial peripheral PTA size matrix that facilitates broader treatment of lesions and delivers exceptional performance.

Ordering Information (Packaging: 1 unit/box)



BALLOON DIAMETER (mm)	BALLOON LENGTH (mm)	RATED BURST PRESSURE (atm)	SHEATH COMPATIBILITY (Fr)	SHAFT LENGTH (190cm)	ORDER QTY
4	100	18	4	R51004010L	
4	150	14	4	R51004015L	
4	200	14	4	R51004020L	
4	250	14	4	R51004025L	
4	300	14	4	R51004030L	
5	20	16	4	R51005002L	
5	30	16	4	R51005003L	
5	40	16	4	R51005004L	
5	60	16	4	R51005006L	
5	80	16	4	R51005008L	
5	100	16	4	R51005010L	
5	150	14	4	R51005015L	
5	200	14	4	R51005020L	
5	250	14	4	R51005025L	
5	300	14	4	R51005030L	
6	20	16	4	R51006002L	
6	30	16	4	R51006003L	
6	40	16	4	R51006004L	
6	60	16	4	R51006006L	
6	80	16	5	R51006008L	
6	100	16	5	R51006010L	
6	150	14	5	R51006015L	
6	200	14	5	R51006020L	
6	250	14	6	R51006025L	
6	300	14	6	R51006030L	
7	20	14	4	R51007002L	
7	30	14	4	R51007003L	
7	40	14	5	R51007004L	
7	60	14	5	R51007006L	
7	80	14	5	R51007008L	
7	100	14	5	R51007010L	
8	20	12	4	R51008002L	
8	30	12	4	R51008003L	
8	40	12	4	R51008004L	
8	60	12	5	R51008006L	
8	80	12	5	R51008008L	
8	100	12	5	R51008010L	
9	20	10	5	R51009002L	
9	30	10	5	R51009003L	
9	40	10	5	R51009004L	
9	60	10	5	R51009006L	
9	80	10	5	R51009008L	
9	100	10	5	R510090010L	
10	20	10	5	R51010002L	
10	30	10	5	R51010003L	
10	40	10	5	R51010004L	
10	60	10	5	R51010006L	
10	80	10	6	R51010008L	
10	100	10	6	R510100010L	

## PRODUCT ORDER FORM

### ZEPHYR® Vascular Compression Band

Designed to provide security, comfort, and simplicity.

Ordering Information (Packaging: 5 units/box)



SIZE	LENGTH (cm)	SKU	ORDER QTY
Large	30	190102	

### RAIN SHEATH® Transradial Thin-Walled Sheath

Next generation thin-walled sheath designed for consistent performance in radial cases.

Ordering Information (Packaging: 5 units/box)



SHEATH SIZE	CANNULA LENGTH	MINI-WIRE	MINI-WIRE COMPATIBILITY	NEEDLE	SKU	ORDER QTY
4F	10cm	POLYMER ANGLED-TIP	0.021" X 45cm	20G IV CATHETER NEEDLE	506410P20G	
		BARE STAINLESS STEEL	0.021" X 45cm	21G BARE NEEDLE	506410S	
		POLYMER ANGLED-TIP	0.021" X 45cm	22G IV CATHETER NEEDLE	506410P22G	
		BARE NITINOL	0.021" X 43cm	21G BARE NEEDLE	506410N	
	16cm	POLYMER ANGLED-TIP	0.021" X 45cm	20G IV CATHETER NEEDLE	506416P20G	
		BARE STAINLESS STEEL	0.021" X 45cm	21G BARE NEEDLE	506416S	
		POLYMER ANGLED-TIP	0.021" X 45cm	22G IV CATHETER NEEDLE	506416P22G	
		BARE NITINOL	0.021" X 43cm	21G BARE NEEDLE	506416N	
5F	10cm	POLYMER ANGLED-TIP	0.021" X 45cm	20G IV CATHETER NEEDLE	506510P20G	
		BARE STAINLESS STEEL	0.021" X 45cm	21G BARE NEEDLE	506510S	
		POLYMER ANGLED-TIP	0.021" X 45cm	22G IV CATHETER NEEDLE	506510P22G	
		BARE NITINOL	0.021" X 43cm	21G BARE NEEDLE	506510N	
	16cm	POLYMER ANGLED-TIP	0.021" X 45cm	20G IV CATHETER NEEDLE	506516P20G	
		BARE STAINLESS STEEL	0.021" X 45cm	21G BARE NEEDLE	506516S	
		POLYMER ANGLED-TIP	0.021" X 45cm	22G IV CATHETER NEEDLE	506516P22G	
		BARE NITINOL	0.021" X 43cm	21G BARE NEEDLE	506516N	
6F	10cm	POLYMER ANGLED-TIP	0.021" X 45cm	20G IV CATHETER NEEDLE	506610P20G	
		BARE STAINLESS STEEL	0.021" X 45cm	21G BARE NEEDLE	506610S	
		POLYMER ANGLED-TIP	0.021" X 45cm	22G IV CATHETER NEEDLE	506610P22G	
		BARE NITINOL	0.021" X 43cm	21G BARE NEEDLE	506610N	
	16cm	POLYMER ANGLED-TIP	0.021" X 45cm	20G IV CATHETER NEEDLE	506616P20G	
		BARE STAINLESS STEEL	0.021" X 45cm	21G BARE NEEDLE	506616S	
		POLYMER ANGLED-TIP	0.021" X 45cm	22G IV CATHETER NEEDLE	506616P22G	
		BARE NITINOL	0.021" X 43cm	21G BARE NEEDLE	506616N	
7F	10cm	POLYMER ANGLED-TIP	0.021" X 45cm	20G IV CATHETER NEEDLE	506710P20G	
		BARE STAINLESS STEEL	0.021" X 45cm	21G BARE NEEDLE	506710S	
		POLYMER ANGLED-TIP	0.021" X 45cm	22G IV CATHETER NEEDLE	506710P22G	
		BARE NITINOL	0.021" X 43cm	21G BARE NEEDLE	506710N	
	16cm	POLYMER ANGLED-TIP	0.021" X 45cm	20G IV CATHETER NEEDLE	506716P20G	
		BARE STAINLESS STEEL	0.021" X 45cm	21G BARE NEEDLE	506716S	
		POLYMER ANGLED-TIP	0.021" X 45cm	22G IV CATHETER NEEDLE	506716P22G	
		BARE NITINOL	0.021" X 43cm	21G BARE NEEDLE	506716N	

## PRODUCT ORDER FORM

### STORQ® STEERABLE GUIDEWIRE

Ordering Information (Packaging: 5 units/box)

PRODUCT	SIZE (in)	LENGTH (cm)	SHAPE	SKU	ORDER QTY
STORQ®	0.035	300	Straight	503456X	
STORQ®	0.035	300	Angled	503456Y	

### SV STEERABLE GUIDEWIRE

Ordering Information (Packaging: 5 units/box)

PRODUCT	SIZE (in)	LENGTH (cm)	SHAPE	SKU	ORDER QTY
SV5	0.018	300	Straight	503558X	
SV8	0.018	300	Straight	503658X	

### STABILIZER® STEERABLE GUIDEWIRE

Ordering Information (Packaging: 5 units/box)

PRODUCT	SIZE (in)	LENGTH (cm)	SHAPE	SKU	ORDER QTY
Support Guidewire	0.014	300	Straight	507300S	
Extra Support Guidewire	0.014	300	Straight	527300E	

### INFINITI® DIAGNOSTIC CATHETER

Non Coated Diagnostic Catheters

Ordering Information (Packaging: 5 units/box)

LENGTH (cm)	SHEATH FIT (F)	SHAPE	SKU	ORDER QTY
110	4	Pig	538450S	
110	5	Pig	534550S	
110	6	Pig	534650S	

### TEMPO AQUA® CATHETER

Hydrophilic Coated Diagnostic Catheters

Ordering Information (Packaging: 1 unit/box)

LENGTH (cm)	SHEATH FIT (F)	SHAPE	SKU	ORDER QTY
125	4	Straight	SRD6667	
125	4	H1	SRD6107	
125	4	Vert	SRD6108	
125	5	Straight	SRD6669	
125	5	H1	SRD6721	
125	5	SIM	SRD6109	