

SLEEK® OTW


.014" Over-the-Wire PTA Catheter

Excellent Pushability and a Small Crossing Profile.


- Choose from 1.25-5 mm diameters and 15-280 mm lengths
- Long 150 cm working catheter lengths
- All balloon sizes compatible with 4F sheaths

PRODUCT BENEFITS:


- **Excellent pushability** with a **small crossing profile** to facilitate the treatment of the most challenging and distal infrapopliteal lesions¹
- Offers a **broad portfolio** enabling you to treat patients with severe vascular disease with fewer inflations



Excellent Pushability¹ (Higher is better)



Small Crossing Profile¹ (Lower is better)



SLEEK[®] OTW


.014" PTA Dilatation Catheter

150 cm Delivery System

ORDERING INFORMATION

Catalog Number	Balloon Diameter (mm)	Balloon Length (mm)	Rated Burst Pressure (atm)	Sheath Fit (F)
4261201X	1.25	15	14	4
4261501X	1.5	15	16	4
4261502X	1.5	20	16	4
4261504X	1.5	40	16	4
4261508X	1.5	80	16	4
4261510X	1.5	100	16	4
4261512X	1.5	120	16	4
4262002X	2	20	16	4
4262004X	2	40	16	4
4262008X	2	80	15	4
4262010X	2	100	15	4
4262012X	2	120	15	4
4262015X	2	150	15	4
4262022X	2	220	15	4
4262028X	2	280	15	4
4262502X	2.5	20	16	4
4262504X	2.5	40	16	4
4262508X	2.5	80	15	4
4262510X	2.5	100	15	4
4262512X	2.5	120	15	4
4262515X	2.5	150	15	4
4262522X	2.5	220	15	4
4262528X	2.5	280	15	4
4263002X	3	20	16	4
4263004X	3	40	16	4
4263008X	3	80	15	4
4263010X	3	100	15	4
4263012X	3	120	15	4

Catalog Number	Balloon Diameter (mm)	Balloon Length (mm)	Rated Burst Pressure (atm)	Sheath Fit (F)
4263015X	3	150	15	4
4263022X	3	220	15	4
4263028X	3	280	15	4
4263502X	3.5	20	16	4
4263504X	3.5	40	16	4
4263508X	3.5	80	15	4
4263510X	3.5	100	15	4
4263512X	3.5	120	15	4
4263515X	3.5	150	15	4
4263522X	3.5	220	15	4
4263528X	3.5	280	15	4
4264002X	4	20	16	4
4264004X	4	40	16	4
4264008X	4	80	15	4
4264010X	4	100	15	4
4264012X	4	120	15	4
4264015X	4	150	15	4
4264022X	4	220	15	4
4264028X	4	280	15	4
4265002X	5	20	14	4
4265004X	5	40	14	4
4265008X	5	80	13	4
4265010X	5	100	13	4
4265012X	5	120	10	4
4265015X	5	150	10	4
4265022X	5	220	10	4
4265028X	5	280	10	4



*There are two marker bands denoting the ends of the working surface on all sizes except the 1.25 mm by 15 mm and 1.5 mm by 15 mm, where a single marker band is centered on the working surface.

CAUTION: Federal (USA) law restricts this device to sale by or on the order of a physician.

Important information: Prior to use, refer to the instruction for use supplied with this device for indications, contraindications, side effects, suggested procedure, warnings and precautions.

CORDIS, Cordis LOGO, and SLEEK are trademarks of Cordis and may be registered in the US and/or in other countries.

All other marks are the property of their respective owners.

© 2021 Cordis. All Rights Reserved. 100591836 09/2021

Cordis